

BELVIR Workshop on IFN

PARTICIPANTS

Seminar

> 20 members of BELVIR labs attended the seminar on Thursday morning
+ extra persons

Workshop

Tim Dierckx	KUL Debyser
Olivier Vandersleyen	Coda-Cerva
Fiona Ingraio	Coda-Cerva
Thu Giang Nguyen	Coda-Cerva
Jochen Lamote	U-Ghent
Niraj Mishra	KUL Neyts
Joanna Zmurko	KUL Neyts
Yuan Gao	ULG Vanderplasschen
Lea Morvan	ULG Vanderplasschen
(Nicolas Gillet) Excused	ULG Willems

Contributors from host lab UCL Michiels:

Michael Peeters
Melissa Drappier
Aurélie De Cock
Muriel Minet
Stephane Messe

IFN workshop



February 12-13th, 2014

All welcome

Auditorium de Visscher,
de Duve Institute,
74 avenue Hippocrate, B-1200 Brussels

9h00-9h30: Registration

9h30-9h55: Welcome session and Introduction on IFN - Thomas Michiels

9h55-10h20: RNase L activity and its inhibition by Theiler's virus - Melissa Drappier

10h20-10h40: Activation and inhibition of PKR - Fabian Borghese

10h40-11h00: Coffee Break

11h00-11h25: HIV triggers a cGAS-dependent type 1 IFN response in primary CD4⁺ T cells that is regulated by Vpu and Vpr - Jolien Vermeire* (U Ghent)

* kindly replaced by B. Verhasselt for medical reasons.

11h20-11h45: Methods to measure IFN - Thomas Michiels

12h00-13h00: Seminar

Fine control of type I interferons signaling : relevance in health and disease

Sandra Pellegrini, Pasteur Institute, Paris

13h00-14h30: Lunch

AFTERNOON and FRIDAY: PRACTICAL TRAINING

Mx immunolabeling


Thursday :

- deparaffination : successive bath of Xylene, Ethanol -> H₂O to rehydrate the sections (5 min each in : Xylene, Xylene again, Ethanol 100, Ethanol 95, Ethanol 70, Ethanol 50, H₂O)
- epitope retrieval : 30 min. 95°C in citrate buffer pH 4.5
- permeabilisation : 5 min in 0,2 % Triton X-100
- blocking : incubation in 100 ul of blocking buffer (goat serum 1/50)
- Primary antibody : diluted 1/150 in blocking solution (overnight 4°C) – wet room(box)

Friday :

- wash : 3 x 5 min in PBS – 0,1% Tween-20
- incubate 1h with secondary antibody : Alexa488-conjugated anti-rabbit antibody (1/400).
- (optional) add Hoechst nuclear stain 1 ug/ml
- wash : 3 x 5 min in PBS – 0,1% Tween-20

Mount with a hydrophilic mounting medium (Mowiol).
Observe under the fluorescent microscope.



IFN-α induces Mx1 expression in the lamina propria
IFN-λ induces Mx1 expression in the epithelium.

IFN assay

Tuesday

- Cells : BALB/3T3 (mouse fibroblasts) and 2ftgh (human fibrosarcoma) cells
- seed 8000 cells per well in 150 ul of medium (DMEM-10% FCS) of in 96-well plates

Wednesday


- Add 150 ul of medium containing IFN (or control) in the first raw of wells.
- make 2-fold dilution series by pipetting 150 ul from well to well (change tips)
- control cells treated with culture medium
- positive control with reference IFN

Thursday

- infect wells with 0,5 PFU per cell of Mengovirus (EMCV) or VSV-GFP
- control non-infected wells

Friday

- observe cytopathic effect under the light microscope
- observe VSV-GFP-infected cells under the fluorescence microscope.
- fix the cells by adding 40 ul formaldehyde 37% per well
- after > 2h, discard the supernatant and color the plates with violet cristal.



RNase L activity

Wednesday

- co-transfection 6-cm dishes (27 cm²) of 293T cells (60-70% confluence) with plasmids expressing
 - muRNase L
 - muRNase L and L*

Thursday

- wash cells with ice-cold PBS
- lyse cells
- scrape cells and centrifuge the lysate
- divide the lysate supernatant into aliquots
- add 2-5A
- Prepare RNA after 5 or 20 min incubation

RNA prep : (according to Chomczynski and Sacchi, 1987)

- add 400 ul of D solution (Guanidine Thiocyanate 4M, Na citrate 25mM, pH7 ; N-lauroylsarcosine 0.5% ; β -mercaptoethanol 0.1M)
- mix vigorously
- add 500 ul of phenol+Na acetate and mix vigorously
- add 125 ul of CHCl₃ and mix vigorously
- centrifuge for 20 min at 4°C
- collect 400 ul supernatant and add 400 ul isopropanol -> Mix -> incubate ON at -20°C

Friday

- centrifuge RNAs for 20 min at 4°C
- wash pellets with icecold Ethanol
- resuspend RNAs in RNase-free water
- run on standard 1 % agarose gels

Observation : unfortunately, degradation was hardly visible in this experiment.

RT-qPCR for IFN and ISG mRNA expression

Miqe guidelines about design and report of qPCR reactions

<http://miqe.gene-quantification.info/>

Web site for primer design

<http://primer3plus.com/cgi-bin/dev/primer3plus.cgi>

Several other web sites exist to define primers and probes.

Note that Taqman probes may be chosen with a higher melting temp (5 to 10°C higher than primers T_m)

Thursday

- mix the qPCR master mix (Sybr Green) with primers

- Primers for IFN-beta (mouse)

Fwd : ATG AAC AAC AGG TGG ATC CTC C
Rev : AGG AGC TCC TGA CAT TTC CGA A

- Primers for Oasl2 (mouse)

Fwd : GGA TGC CTG GGA GAG AAT CG
Rev : TCG CCT GCT CTT CGA AAC TG

- prepare the plasmid dilutions for the standard curve

- prepare the RT-qPCR plate :

- 5 ul of diluted cDNA

- 20 ul of qPCR Mix

- controls include :

- non-infected samples

- cDNA samples generated without RT

- H₂O

- Run (note that more recent enzymes allow to do much faster runs)

95°C – 10 min (to activate « hot-start » Taq polymerase)

-> 40 cycles :

95°C – 15 sec : denaturation

60°C – 1 min : primer annealing – elongation – fluorescence measurements

-> 30 x 10 sec steps of + 0,5°C : melting curve

(melting curve allows to assess the specificity of the products
when Sybr green is used instead of Taqman probes)

Friday

- analysis of the data